

SURGICAL LOUPES

BY SEILER MEDICAL DIVISION



2.5x, 3.0x and 3.5x magnifications. Provides vision not attainable with the naked eye.

HIGH RESOLUTION – highest grade glass for excellent clarity.

WATER SEALED – optics are sealed to prevent dust and moisture.

DOUBLE HINGE FLIP-UP DESIGN – for superior optical alignment.

GREAT DEPTH OF FIELD – comfortably extends the natural accommodation capability of the human eye.

WIDE-FIELD – highly corrected optics that allows the entire mouth to be viewed.

WORKING DISTANCE – Short (340mm or 13.4 inches), Regular (420mm or 16.5 inches), Long (500mm or 19.7 inches).

All Seiler loupes allow the user to set your interpupillary, convergence and declination angle for perfect optical alignment.

FRAME STYLES



Lightweight titanium frames allow the user to adapt an RX lens to the loupes.

Available Colors: ■ Metallic Blue ■ Metallic Silver
■ Metallic Red ■ Black



The Sport Frame is a large durable frame that offers added protection against splash-back and debris.

Available Colors: ■ Metallic Silver ■ Black

MAGNIFICATION LEVELS

2.5x

The entry level magnification that provides a wide field of view, which is exceptionally easy to get used to and is the most common magnification bought.

3.0x

The intermediate magnification level of which is good for a user that is familiar with loupes or has a need for additional magnification not provided by the 2.5x power.

3.5x

The most powerful loupe level in the Seiler Line of Loupes. This magnification level is used for specialists and users that have a smaller working field.

WORKING DISTANCE

short | 13.4" (340mm)

The short working distance is 13.4 inches or 340mm.

regular | 17" (420mm)

The regular working distance is 16.5 inches or 420mm. The regular working distance is the most common working distance.

Long | 19" (500mm)

The long working distance is 19.7 inches or 500mm. The long working distance is generally used by a tall person.

LED HEADLIGHT

The Seiler LED Loupe Light is designed with precision engineering to provide the brightest light available. Additional features are:

- Ultra Lightweight Design.
- Over 50,000 LUX.
- Supplies bright, white, shadow-free light.
- Variable Illumination
- Rechargeable Battery Pack makes this a complete portable, compact unit that allows the user freedom of movement.
- Seiler's high quality LED Light operating life expectancy is approximately 500 charge/discharge cycles.



Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
+61 2 8090 0994
info@mintdevices.com.au
www.mintdevices.com.au

**SEILER**
MEDICAL DIVISION
info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au