




PROMISEVISION 3D™ SURGICAL MICROSCOPE

Dental



3D SURGICAL MICROSCOPE TECHNICAL SPECIFICATIONS

MOUNTING OPTIONS

- Floor
- Fixed Floor
- Ceiling
- Wall
- High Wall
- Tabletop

- 4.5X - 28X Magnification Range
- Integrated Movement of Cameras and Monitor for Freedom of Movement Around the Patient
- Monitor and Cameras are Linked to Maintain the Most Understandable Frame of Reference, Properly Oriented for the Dentist to Work
- Seamless, Magnified Video Image at 60 FPS
- Fastest Image Processing Available - No Latency in the Video Image
- Full HD Monitor for Best Possible 3D Image Resolution
- Unprecedented Depth of Field (up to 60 mm) and Field of View (Up to 120mm for Full Arch View)
- "Head Up" 3D View Provides Ergonomic Posture (Healthy Posture) for Comfort and Productivity
- Natural Intuitive 3D Image with Normal Depth Perception and No Loss of Clarity
- Elimination of the Binocular Eyepieces Improves Productivity and Patient Awareness
- Coaxial LED Illumination (5,000°K)



Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
+61 2 8090 0994
info@mintdevices.com.au
www.mintdevices.com.au

 **SEILER**
MEDICAL DIVISION
info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au