



# ALPHA AIR SERIES

FOCUS ON PRECISION.



# DENTAL OPERATING MICROSCOPE



## ALPHA AIR SERIES

- 6 step turret magnification system (3 step turret available)
- LED Illumination + 150,000 LUX, 50,000 hour life
- Widefield eyepieces with dioptic locks: 10x, 12.5x 16x, and 20x
- True Coaxial Illumination="Shadow Free"
- Objective Lens: 175mm, 200mm, 250mm, 300mm, and 400mm  
OPTIONAL: Vario Focus Lens - 200mm-350mm
- 0-220° Inclinal Binocular Head (Fixed 45° head available)
- Composite and Red Free Filter standard
- USB and 12V DC power support directly on optic pod (charge Handycams, GoPro, DSLR, HD CCD)

## MOUNTING OPTIONS:


Floor | Ceiling | Wall | High Wall | Tabletop



## ACCESSORIES:

- 50/50-50/50, 50/50-20/80 Beamsplitter
- Full 1080P HD Live Video
- Assistant Binocular Head
- Extender/Rotation Device
- Laser Filters
- Iris Diaphragm
- Vario Focus
- All DSLR, Handycam, and GoPro Camera attachments

**MINTDEVICES**



"Leaders in Minimally Invasive Technology"

Mint Devices (Australia) Pty Ltd  
PO Box 290 Potts Point NSW 2011  
+61 2 8090 0994  
info@mintdevices.com.au  
www.mintdevices.com.au



info@seilermicroscopes.com.au  
www.seilermicroscopes.com.au  
www.mintdevices.com.au