


MICROLUX IV



MINTDEVICES



"Leaders in Minimally Invasive Technology"

 **SEILER**
PRECISION MICROSCOPES

MICROLUX IV



Microlux IV (MLX817)


Plan Achromat Binocular Microscope

The MicroluxIV is our most versatile compound microscope with optional video or digital configurations as well as phase contrast, LED, and MOHS for dermatology procedures

Specifications:

- 5 Infinity Corrected Plan Objectives: 4x, 10x, 20x, 40x, 100x
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Kohler Illumination System: Halogen 12V-20W (LED optional)
- Eyepiece: WF10X/20mm
- interpupillary Distance: 50-75mm
- Coaxial Coarse and Fine Focusing Control
- Nosepiece: Reversed Quintuple Click-Stop Revolving Mechanism with Multiple Ball Bearings
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Double Layer Mechanical Stage with Removable Slide Holder, 190x140mm
- Travel Stage Area 78x55mm with A Right-Hand Stage Handle
- Anti-Fungal Systems
- 110V-240V Universal Power Supply
- Limited Lifetime Warranty

MINTDEVICES



"Leaders in Minimally Invasive Technology"

Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
+61 2 8090 0994
info@mintdevices.com.au
www.mintdevices.com.au

 **SEILER**
MEDICAL DIVISION

info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au