

Microlux IV for

M

O

H

S

Dermatology

SPECIFICATIONS:

- Infinity Plan Objectives: 2.5x, 4x, 10x, 40x, 100x
- Kohler Halogen Illumination system, 6V/20W
- Eyepiece: WF10X/20
- Coaxial Fine Focusing adjustment with rack and pinion mechanism
- Nosepiece: Quintuple click-stop revolving mechanism with multiple ball bearings
- Anti-fungal systems
- Seidentopf Type Binocular head, 30° inclined, 360° rotatable
- Double layer mechanical stage with removable slide holder, 190x140mm




"This has been a great microscope for us. There are some colors and slides that are just easier to read on my Seiler microscope. It works very well for Mohs. This used to be a backup microscope and has become my go to microscope."

*~ Dr. Stephen Craig, MD
North Idaho Dermatology*

 **SEILER**
PRECISION MICROSCOPES

MINTDEVICES



"Leaders in Minimally Invasive Technology"



Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
info@mintdevices.com.au
www.mintdevices.com.au
+61 2 8090 0994

info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au

 **SEILER**
PRECISION MICROSCOPES