


# SEILERSCOPE



MINTDEVICES



"Leaders in Minimally Invasive Technology"

 **SEILER**  
PRECISION MICROSCOPES

# SEILERSCOPE




## SeilerScope (SXS820)

Achromat Binocular Microscope

### Specifications:

- Achromat Objectives: 4, 10, 40, 100x (oil)
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepieces: WF10X/18mm
- 6V20W Halogen Illumination
- Coaxial Coarse And Calibrated Fine Focus Control
- Graduated, Right Hand Double Plated Mechanical Stage
- Reversed Four-Hole Nosepiece
- Focusing Tension Adjustment Knob
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- 110V or 220V switchable power supply
- Anti-Fungal System
- Optional Accessories: Live Video, Plan Objectives
- Limited Lifetime Warranty

MINTDEVICES



*"Leaders in Minimally Invasive Technology"*

Mint Devices (Australia) Pty Ltd  
PO Box 290 Potts Point NSW 2011  
+61 2 8090 0994  
[info@mintdevices.com.au](mailto:info@mintdevices.com.au)  
[www.mintdevices.com.au](http://www.mintdevices.com.au)

 **SEILER**

PRECISION MICROSCOPES

[info@seilmicroscopes.com.au](mailto:info@seilmicroscopes.com.au)  
[www.seilmicroscopes.com.au](http://www.seilmicroscopes.com.au)  
[www.mintdevices.com.au](http://www.mintdevices.com.au)