

# WESTLAB III



MINTDEVICES



"Leaders in Minimally Invasive Technology"

 **SEILER**

PRECISION MICROSCOPES

# WESTLAB III




## Westlab III (WLB814)

Semi-Plan Binocular Microscope

Includes:

- Semi-Plan Objectives: 4, 10, 40, 100X oil
- 6V 20W Halogen Illumination
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepiece: WF10X/20mm
- Interpupillary distance: 50-75mm
- Reversed Quadruple click-stop nosepiece
- Coaxial Coarse And Calibrated Fine Focus Control
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Mechanical stage with removable slide holder 142 x 140mm; travel area with a right-hand stage handle
- Anti-Fungal System
- 100V – 240V universal power supply
- Optional: Live Video Feature
- Limited Lifetime Warranty

## MINTDEVICES



*"Leaders in Minimally Invasive Technology"*

Mint Devices (Australia) Pty Ltd  
PO Box 290 Potts Point NSW 2011  
+61 2 8090 0994  
[info@mintdevices.com.au](mailto:info@mintdevices.com.au)  
[www.mintdevices.com.au](http://www.mintdevices.com.au)

 **SEILER**  
PRECISION MICROSCOPES  
[info@seilmicroscopes.com.au](mailto:info@seilmicroscopes.com.au)  
[www.seilmicroscopes.com.au](http://www.seilmicroscopes.com.au)  
[www.mintdevices.com.au](http://www.mintdevices.com.au)