



ALPHA AIR SERIES

FOCUS ON PRECISION.



ENT OPERATING MICROSCOPE



ALPHA AIR SERIES

- 6 step turret magnification system (3 step turret available)
- LED Illumination + 150,000 LUX, 50,000 hour life
- Widefield eyepieces with dioptic locks: 10x, 12.5x 16x, and 20x
- True Coaxial Illumination="Shadow Free"
- Objective Lens: 175mm, 200mm, 250mm, 300mm, and 400mm
OPTIONAL: Vario Focus Lens - 200mm-350mm
- Straight Head Standard (0-220° Inclunable Binocular Head or Inclined Fixed 45° head available)
- Red Free Filter standard
- USB and 12V DC power support directly on optic pod (charge Handycams, GoPro, DSLR, HD CCD)

MOUNTING OPTIONS:


Floor | Ceiling | Wall | High Wall | Tabletop



ACCESSORIES:

- 50/50-50/50, 50/50-20/80 Beamsplitter
- Assistant Binocular Head
- Extender/Rotation Device
- Laser Filters
- Iris Diaphragm
- Vario Focus
- DSLR Camera (Nikon, Canon, or Sony), Full 1080P HD Live Video, and GoPro

MINTDEVICES



"Leaders in Minimally Invasive Technology"

Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
+61 2 8090 0994
info@mintdevices.com.au
www.mintdevices.com.au

SEILER
MEDICAL DIVISION

info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au