

MINTDEVICES



"Leaders in Minimally Invasive Technology"




EVOLUTION ZOOM

FOCUS ON PRECISION



MULTIDISCIPLINARY SURGICAL MICROSCOPE



Neurosurgery, Spine, ENT, and Plastic



TECHNICAL SPECIFICATIONS:


- Magnification Changer: ZOOM 1:6 Motorized
- APOCHROMATIC OBJECTIVES f=175mm-400mm
- Widefield Eyepieces 10x, 12.5x, and 20x, Diopter adjustment +6/-6
- LED Illumination + 100,000 LUX, 50,000 hour life
OPTIONAL: Xenon Illumination=180,000 LUX
Halogen Illumination=80,000 LUX
- Motorized 40mm Fine Focus
- Optional XY Unit range 50mm with automatic reset Button
- Control Unit: Multi-functional Display for illumination, focus, zoom, and automatic reset
- Hand/Foot Switch: 14 function foot pedal and 8 function hand control grip

PN: EVO-100Z
PN: EVO-100XY

ACCESSORIES:

- 50/50 -50/50, 50/50-20/80 Beamsplitter
- Full 1080P HD Live Video
- All DSLR Camera Attachments
- Assistant Binocular Teaching Head
- Extender/Rotation Device
- Laser Filters
- Iris Diaphragm

MINTDEVICES



"Leaders in Minimally Invasive Technology"

Mint Devices (Australia) Pty Ltd
PO Box 290 Potts Point NSW 2011
+61 2 8090 0994
info@mintdevices.com.au
www.mintdevices.com.au

 **SEILER**
MEDICAL DIVISION

info@seilermicroscopes.com.au
www.seilermicroscopes.com.au
www.mintdevices.com.au